

AD-A114 811

LETTERMAN ARMY INST OF RESEARCH PRESIDIO OF SAN FRANC--ETC F/6 6/5  
COMPUTER ASSISTED LABELLING IN MUTAGENICITY TESTING. I. THE DRO--ETC(U)  
MAR 82 W W JEDERBERG, R A WIRTZ, N R POWERS  
LAIR-82-32TN

UNCLASSIFIED

NL

1 OF 1  
50  
3750



END  
DATE  
FILMED  
6 82  
DTIC

C12

AD A114511

DTIC FILE COPY

| REPORT DOCUMENTATION PAGE   |                       | READ INSTRUCTIONS<br>BEFORE COMPLETING FORM   |
|---|-----------------------|---|
| 1. REPORT NUMBER  | 2. GOVT ACCESSION NO. | 3. RECIPIENT'S CATALOG NUMBER   |
|   | AD-A114511            |   |
| 4. TITLE (and Subtitle)<br>Computer Assisted Labelling In Mutagenicity<br>Testing: I. The <u>Drosophila Melanogaster</u><br>Sex-Linked Recessive Lethal Assay.  |                       | 5. TYPE OF REPORT & PERIOD COVERED<br>Final Aug-80 to Aug-81                                    |
|   |                       | 6. PERFORMING ORG. REPORT NUMBER  |
| 7. AUTHOR(s)<br>Warren W. Jederberg, MS, CPT, MS<br>Robert A. Wirtz, PHD, CPT, MS<br>Nelson R. Powers, PHD, CPT, MS   |                       | 8. CONTRACT OR GRANT NUMBER(s)  |
| 9. PERFORMING ORGANIZATION NAME AND ADDRESS<br>Toxicology Services Group<br>Division of Research Support<br>Letterman Army Institute of Research<br>Presidio of San Francisco, CA 94129   |                       | 10. PROGRAM ELEMENT, PROJECT, TASK<br>AREA & WORK UNIT NUMBERS<br>Project 35162772.875<br>WU302 |
| 11. CONTROLLING OFFICE NAME AND ADDRESS<br>US Army Medical Research and Development Command<br>Fort Detrick, Frederick, MD 21701  |                       | 12. REPORT DATE<br>March 1982*  |
|   |                       | 13. NUMBER OF PAGES<br>28   |
| 14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)<br>US Army Medical Research and Development Command<br>Fort Detrick, Frederick, MD 21701  |                       | 15. SECURITY CLASS. (of this report)<br>UNCLASSIFIED  |
|   |                       | 15a. DECLASSIFICATION/DOWNGRADING<br>SCHEDULE   |
| 16. DISTRIBUTION STATEMENT (of this Report)<br>This document has been approved for public release and sale, its distribution<br>is unlimited.   |                       |   |
| 17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)  |                       |   |
| 18. SUPPLEMENTARY NOTES<br>*LAIR Technical Note No. 82-32TN   |                       |   |
| 19. KEY WORDS (Continue on reverse side if necessary and identify by block number)<br>Drosophila melanogaster<br>Sex-linked recessive lethal assay<br>Mutagenicity testing  |                       |   |
| 20. ABSTRACT (Continue on reverse side if necessary and identify by block number)<br>This report presents the FORTRAN V programs utilized for the production of labels<br>and cards for the <u>Drosophila melanogaster</u> sex-linked recessive lethal mutageni-<br>city assay. This program provides for the efficient and correct labelling of<br>articles used in the assay. The use of this program to produce parallel labels<br>and cards not only reduces preparation time, but assures the proper sequence<br>and completeness of the assay. This along with data storage and archiving also<br>insures compliance with the requirements of the Food and Drug Administrations<br>Good Laboratory Practices Act. |                       |   |

DTIC  
ELECTE  
MAY 18 1982  
H

*TECHNICAL NOTE NO. 82-32TN*

**COMPUTER-ASSISTED LABELLING IN MUTAGENICITY TESTING**

**I. The *Drosophila Melanogaster* Sex-Linked Recessive Lethal Assay**

*WARREN W. JEDERBERG, MS, CPT MS*

*ROBERT A. WIRTZ, PhD, CPT MS*

*and*

*NELSON R. POWERS, PhD, CPT MS*

**DIVISION OF CUTANEOUS HAZARDS**

**TOXICOLOGY GROUP,**

**DIVISION OF RESEARCH SUPPORT**

**MARCH 1982**

**LETTERMAN ARMY INSTITUTE OF RESEARCH  
PRESIDIO OF SAN FRANCISCO, CALIFORNIA 94129**

COMPUTER-ASSISTED LABELLING IN MUTAGENICITY TESTING.

I. The *Drosophila Melanogaster* Sex-Linked Recessive Lethal Assay

Reproduction of this document in whole or in part is prohibited except with the permission of the Commander, Letterman Army Institute of Research, Presidio of San Francisco, California 94129. However, the Defense Technical Information Center is authorized to reproduce the document for United States Government purposes.

Destroy this report when it is no longer needed. Do not return it to the originator.

Citation of trade names in this report does not constitute an official endorsement or approval of the use of such items.

This material has been reviewed by Letterman Army Institute of Research and there is no objection to its presentation and/or publication. The opinions or assertions contained herein are the private views of the author(s) and are not to be construed as official or as reflecting the views of the Department of the Army or the Department of Defense. (AR 360-5)

*John H. Marshall* 10 March 1962  
(Signature and date)

# PREFACE

This is the first in a series of technical notes on the utilization of the computer facilities at Letterman Army Institute of Research to assist in mutagenicity testing as part of the Institute's toxicology program. These reports will detail the use of the computer in labelling articles, recording, storing and retrieving data, and data analysis. Both the Salmonella mutagenicity and the Drosophila melanogaster sex-linked recessive lethal mutagenicity assays will be covered.



|                    |                                     |
|--------------------|-------------------------------------|
| Accession For      |                                     |
| NTIS GRA&I         | <input checked="" type="checkbox"/> |
| DTIC TAB           | <input type="checkbox"/>            |
| Unannounced        | <input type="checkbox"/>            |
| Justification      |                                     |
| By                 |                                     |
| Distribution/      |                                     |
| Availability Codes |                                     |
| Dist               | Avail and/or<br>Special             |
| A                  |                                     |

#### ACKNOWLEDGMENTS

The authors express their appreciation to the Information Sciences Group for making computer facilities readily available and to all those who have made suggestions for improving the utility of the program described herein.

## COMPUTER-ASSISTED LABELLING IN MUTAGENICITY TESTING

### I. The *Drosophila Melanogaster* Sex-Linked Recessive Lethal Assay

The mutagenicity testing of materials in compliance with Federal regulations (1) requires the implementation of an extensive labelling system which must be integrated with data collection, storage, data retrieval and analysis. The *Drosophila melanogaster* sex-linked recessive lethal assay is an often used mutagenicity test which requires the unique numbering of the test flies and their progeny, dilutions of test and control materials, and their multiple replicates. The detection of test materials with low mutagenic activity requires the concurrent testing of 7,000-10,000 X-chromosomes, depending on the spontaneous mutation frequency of the stock colonies (2,3), from both negative control and test material treated flies. Due to space and personnel limitations, we currently divide the testing of a single compound into four runs of approximately 2,500 X-chromosomes each for the negative control and test material treated flies as well as 500 X-chromosomes from positive control treated flies (4). This requires the generation of 880 unique, sequential labels and brood cards for testing a single compound. These items are currently hand labelled in most laboratories. We have designed and implemented a FORTRAN V program for the rapid generation of the large numbers of uniquely identified labels and cards. The use of this system has greatly reduced the time spent generating these materials, eliminated errors in numbering, and ensured continuity from the initiation to the termination of the assay.

#### PROGRAM DESCRIPTION

The program (See Appendix A - SLRLABELS.FR) generates the labels and cards for the individual male flies that will be exposed to negative control, positive control, and test compounds.

As seen from the sample run (Figure 1), the program displays the following information (the underlined characters were entered at the time of the program execution by the user): The number of the last run for which labels and cards were prepared and the sequence number of the last male fly for which labels and cards were generated.

As can also be seen in the sample run (Figure 1) and the sample labels and cards (Figures 2-5), the identification number of the males is composed of two fields: the first designates the nature of the compound (C - negative control, P - positive control, and T - test compound) and which of the series of compounds it is by an integer

designator from 1 to 5. In the case of test substances, a repetition number is also printed. The second field of the male identification is the integer sequence number of the fly.

The program requests the following information from the user (refer to Figure 1):

- The initial sequence number to be used for the flies of this run.
- The maximum sequence number to be used for the flies of this run. If the maximum sequence number, minus the minimum sequence number, is greater than 100, then the user must re-enter both values. This feature is utilized because experience has shown that 25 is a reasonable practical limit for setting up the assay in our laboratory.
- The new run number.
- The number of negative control compounds (a maximum of five is allowed) for which labels and cards are to be prepared to permit the concurrent testing of several control compounds.
- The codes for each of the negative control compounds. The user is allowed up to six characters of any type for the designation of each compound.
- The number of positive control compounds (a maximum of five is allowed) for which labels and cards are to be prepared.
- The maximum sequence number for the positive control compounds (usually, there is no need for as many flies to be exposed to the positive control as to the negative control and test compounds).
- The codes for the positive control compounds. The user is allowed up to six characters of any type for the designation of each compound.
- The number of test compounds (a maximum of five is allowed).
- The codes for the test compounds. Again, the user is allowed up to six characters of any type for the designation of each compound.
- The repetition number for each test compound. As is illustrated in Figure 1, no number greater than 99 is allowed.
- The dates for each of the four broods.

As the text for the labels and cards is generated by brood, a message is displayed at the user's terminal to indicate what has been prepared and the files in which they are stored.

After the text for all the labels and cards have been placed in the files DLABELS and BCARDS, the user is so notified and may print the labels and cards out on an appropriate device (usually a printer with a tractor feed). The program is formatted so that the labels and cards should be printed on continuous-feed single-width stock material. These may be obtained from MISCO, 936 Holmdel Keyport Road, Box 399, Holmdel, N.J. 07733, or other distributor of computer related items. The labels are the standard 1 X 3.5 inches and the cards are standard 3 x 5 inches.

All cards and labels are generated so that they are uniquely identified and parallel one another. They are prepared in the sequence in which they will be used. The cards provide permanent records and greatly facilitate data entry into permanent data files on the computer.

#### DISCUSSION

By utilizing the program and subroutines presented in the report we have realized a significant savings in time for the preparation of materials for the *Drosophila melanogaster* sex-linked recessive lethal mutagenicity assay. The labels and cards have also been of great help in maintaining order and uniformity while running this assay.

#### CONCLUSION

None.

#### RECOMMENDATION

None.

#### REFERENCES

1. FOOD AND DRUG ADMINISTRATION. Good laboratory practices regulations. Federal Register 43(163):37336-37403, 1978
2. WURGLER, F.E., U. GRAF and W. BERCHTOLD. Statistical problems connected with the sex-linked recessive lethal test in Drosophila melanogaster. I. The use of the Kastenbaum-Bowman test. Arch Genet 48:158-178, 1975
3. LEE, W.R., S. ABRAHAMSON, R. VALENCIS, E.S. VON HALLE, F.E. WURGLER, and S. ZIMMERLING. The sex-linked recessive lethal test for mutagenesis in Drosophila melanogaster. Mutation Research (In Press)
4. WIRTZ, R.A., N.R. POWERS and J.T. FRUIN. Mutagenicity testing using the Drosophila melanogaster sex-linked recessive lethal assay. Institute Report No. 112, Letterman Army Institute of Research, Presidio of San Francisco, CA 94129, February 1982
5. WIRTZ, R.A. and H.S. Semey. The Drosophila kitchen - equipment, media preparation and supplies. Drosophila Information Services (In Press)

) X SLRLABELS

SLRL-DROSOPHILA LABELS . . . . .

THE LAST RUN WAS: 33  
THE LAST SEQUENCE NUMBER WAS: 813

GLP STUDY NO. ? (15) 12345

INITIAL SEQUENCE NUMBER ? 810  
MAXIMUM SEQUENCE NUMBER ? 813

WHAT IS THE NEW RUN NUMBER ? 33

HOW MANY NEGATIVE CONTROL COMPOUNDS ? (MAX.-5) 1

WHAT IS THE CODE FOR THE NEGATIVE CONTROL #: 1 ? (AN/6)  
NEGONE

HOW MANY POSITIVE CONTROL COMPOUNDS ? (MAX.-5) 1  
MAXIMUM SEQUENCE NUMBER FOR POSITIVE CONTROLS ? 811

WHAT IS THE CODE FOR THE POSITIVE CONTROL #: 1 ? (AN/6)  
POSONE

HOW MANY TEST COMPOUNDS ? (MAX.-5) 2

WHAT IS THE CODE FOR TEST COMPOUND #: 1 ? (AN/6)  
T-ONE

WHAT IS THE REPETITION NUMBER FOR THIS COMPOUND ? (12) 2

WHAT IS THE CODE FOR TEST COMPOUND #: 2 ? (AN/6)  
T-TWO

WHAT IS THE REPETITION NUMBER FOR THIS COMPOUND ? (12) 100  
WHAT IS THE REPETITION NUMBER FOR THIS COMPOUND ? (12) 45

Figure 1. Sample run of SLRLABELS.FR

```

      INPUT BROOD DATES
DATE FOR BROOD: 1 ? (XXMONXX)
99DEC99
DATE FOR BROOD: 2 ? (XXMONXX)
99DEC99
DATE FOR BROOD: 3 ? (XXMONXX)
99DEC99
DATE FOR BROOD: 4 ? (XXMONXX)
99DEC99

```

WAIT FOR OUTPUT . . . . .

RUN: 33

BROOD: 1

```

LABELS AND CARDS FOR C1-NEGONE ARE IN 'DLABELS' AND 'BCARDS'
LABELS AND CARDS FOR P1-POSONE ARE IN 'DLABELS' AND 'BCARDS'
LABELS AND CARDS FOR T1-T-ONE ARE IN 'DLABELS' AND 'BCARDS'
LABELS AND CARDS FOR T2-T-TWO ARE IN 'DLABELS' AND 'BCARDS'

```

BROOD: 2

```

LABELS AND CARDS FOR C1-NEGONE ARE IN 'DLABELS' AND 'BCARDS'
LABELS AND CARDS FOR P1-POSONE ARE IN 'DLABELS' AND 'BCARDS'
LABELS AND CARDS FOR T1-T-ONE ARE IN 'DLABELS' AND 'BCARDS'
LABELS AND CARDS FOR T2-T-TWO ARE IN 'DLABELS' AND 'BCARDS'

```

BROOD: 3

```

LABELS AND CARDS FOR C1-NEGONE ARE IN 'DLABELS' AND 'BCARDS'
LABELS AND CARDS FOR P1-POSONE ARE IN 'DLABELS' AND 'BCARDS'

```

Figure 1. Sample run of SLRLABELS.FR  
(Continued)

LABELS AND CARDS FOR T1-T-ONE ARE IN 'DLABELS' AND 'BCARDS'  
LABELS AND CARDS FOR T2-T-TWO ARE IN 'DLABELS' AND 'BCARDS'

BROOD: 4

LABELS AND CARDS FOR C1-NEGONE ARE IN 'DLABELS' AND 'BCARDS'  
LABELS AND CARDS FOR P1-POSONE ARE IN 'DLABELS' AND 'BCARDS'  
LABELS AND CARDS FOR T1-T-ONE ARE IN 'DLABELS' AND 'BCARDS'  
LABELS AND CARDS FOR T2-T-TWO ARE IN 'DLABELS' AND 'BCARDS'

ALL LABELS AND CARDS ARE READY IN 'DLABELS' AND 'BCARDS'

)

Figure 1. Sample run of SLRLABELS.FR  
(Continued)

GLP STUDY NO. 12345  
RUN: 33  
C1- 813 BR: 1 99DEC99  
COMPOUND CODE: NEGONE  
NOTES:

GLP STUDY NO. 12345  
COMPOUND CODE: NEGONE  
C1- 813 BR: 1 99DEC99 RUN: 33

-----  
F2 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS: \_\_\_\_\_  
FAILED \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----  
F3 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS: \_\_\_\_\_  
FAILED \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----  
NOTES:

Figure 2. Sample Negative Control  
Label and Card

GLP STUDY NO. 12345  
RUN: 33  
P1- 810 BR: 1 99DEC99  
COMPOUND CODE: POSONE  
NOTES:

GLP STUDY NO. 12345  
COMPOUND CODE: POSONE  
P1- 810 BR: 1 99DEC99 RUN: 33

-----  
F2 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS \_\_\_\_\_  
FAILURES \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----  
F3 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS \_\_\_\_\_  
FAILURES \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----  
NOTES:

Figure 3. Sample Positive Control  
Label and Card

GLP STUDY NO. 12345  
RUN: 33 REP. #: 2  
T1- 810 BR: 1 99DEC99  
COMPOUND CODE: T-ONE  
NOTES:

GLP STUDY NO. 12345 REP. #: 2  
COMPOUND CODE: T-ONE  
T1- 810 BR: 1 99DEC99 RUN: 33

-----  
F2 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS: \_\_\_\_\_  
FAILURES \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----  
F3 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS: \_\_\_\_\_  
FAILURES \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----  
NOTES:

Figure 4. Sample Test Compound  
Label and Card

GLP STUDY NO. 12345  
RUN: 33 REP. #: 45  
T2- 810 BR: 1 99DEC99  
COMPOUND CODE: T-TWO  
NOTES:

GLP STUDY NO. 12345 REP. #: 45  
COMPOUND CODE: T-TWO  
T2- 810 BR: 1 99DEC99 RUN: 33

-----  
F2 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS \_\_\_\_\_

FAILURES \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----

F3 CROSS MEDIUM BATCH #: \_\_\_\_\_  
DATE: \_\_\_\_\_ INITIALS \_\_\_\_\_

FAILURES \_\_\_\_\_ LETHALS \_\_\_\_\_ NONLETHALS \_\_\_\_\_  
-----

NOTES:

Figure 5. Sample Test Compound  
Label and Card

|  |    |
|--|----|
| Appendix A. Listing for SRLABELS.FR . . . . .            | 14 |
| Appendix B. Listing for Subroutine NBLABELS.FR . . . . . | 19 |
| Appendix C. Listing for Subroutine NTBLABELS.FR. . . . . | 21 |

APPENDICES

# APPENDIX A

## SLRLABELS.FR

```

C --- PROGRAM NAME:          SLRLABELS.FR
C --- WRITTEN BY:           WARREN W. JEDENBERG
C --- SUBROUTINES:          NBLABELS.FR, NTBLABELS.FR
C --- INPUT:                AT TIME OF EXECUTION:
C ---      RANGE OF LABELS, RUN NUMBER, BROOD DATES
C ---      COMPOUND CODES, COMPOUND NUMBERS
C --- FORMAT FOR INPUT:      CRT OR DASHES
C --- PURPOSE:              TO GENERATE LABELS AND CARLS FOR THE
C ---      SEX-LINKED RECESSIVE LETHAL DROSOPHILA ASSAY.
C ---      THIS PROGRAM AND SUBROUTINES GENERATE THE LABELS FOR
C ---      THE NEGATIVE CONTROL COMPOUND(S), THE POSITIVE CONTROL
C ---      COMPOUND(S), AND THE TEST COMPOUND(S) IN ORDER OF BROOD.
C ---      A MAXIMUM SEQUENCE OF 100 IS ALLOWED AND UP TO 5 OF EACH
C ---      COMPOUND. THE LABELS ARE STORED IN 'DLABELS'.

C --- THE SUBROUTINES NBLABELS.FR AND NTBLABELS.FR ARE USED.
C ---
C --- INTEGER CCOM, PCOM, TCOM, IR, UR, B, TL, IDATE
C --- DIMENSION CCOM(6,5), PCOM(6,5), TCOM(6,5), IDATE(7,4), NREP(5)

C *** LOG ON NOTE
C --- TYPE
C --- TYPE " SERL-DROSOPHILA LABELS . . . . . "
C --- TYPE

C ---
C *** GET THE LAST RUN AND SEQUENCE NUMBERS
C --- OPEN 2, "LDLABEL", ATT="SIB"
C --- READ(2,1) NR,LL
C --- FORMAT (2(1X,I4))
C --- WRITE (10,2) NR,LL
C --- FORMAT(/,1X,"THE LAST RUN WAS:",1X,I4,/,
C --- 1X,"THE LAST SEQUENCE NUMBER WAS:",1X,I4)
C --- CLOSE 2
C ---

C *** GET THE GLP STUDY NUMBER
C --- 100 CONTINUE
C --- TYPE
C --- ACCEPT " GLP STUDY NO. ? (15) ",GSN
C --- IF (GSN.GT.99999) GO TO 100
C ---

C *** INPUT AND CHECK THE SEQUENCE RANGE
C --- 100 CONTINUE
C --- TYPE

```

APPENDIX A  
(CONTINUED)

```

ACCEPT " INITIAL SEQUENCE NUMBER ? ",IR
ACCEPT " MAXIMUM SEQUENCE NUMBER ? ",UR
C --- CHECK RANGE
TL = UR-IR
IF (TL.LE.100) GO TO 300
TYPE
TYPE " ++++++ERROR+++++"
TYPE " MAXIMUM - INITIAL > 100"
GO TO 700
300 CONTINUE
C ---
C *** INPUT NEW RUN NUMBER
TYPE
ACCEPT " WHAT IS THE NEW RUN NUMBER ? ",NR
TYPE
C ---
C *** INPUT COMPOUND CODES
C --- NEGATIVE CONTROL COMPOUND(S)
400 CONTINUE
ACCEPT " HOW MANY NEGATIVE CONTROL COMPOUNDS ? (MAX.-5) ",NNC
IF (NNC.GT.5) GO TO 400
TYPE
DO 90 I=1,NNC
WRITE (10,3) I
3 FORMAT (1X,"WHAT IS THE CODE FOR THE NEGATIVE CONTROL #:",1X,11,1X,
1"? (AN/6)")
READ (11,4) (CCOM(K,I) K=1,6)
4 FORMAT (6A1)
90 CONTINUE
TYPE
TYPE
C ---
C --- POSITIVE CONTROL COMPOUND(S)
500 CONTINUE
ACCEPT " HOW MANY POSITIVE CONTROL COMPOUNDS ? (MAX.-5) ",NPC
IF (NPC.GT.5) GO TO 500
C *** GET NUMBER FOR POSITIVE CONTROL SEQUENCE
900 CONTINUE
ACCEPT " MAXIMUM SEQUENCE NUMBER FOR POSITIVE CONTROLS ? ",PNR
PTL = PNR - IR
IF (PTL.LE.100) GO TO 800
TYPE
TYPE " ++++++ERROR+++++"
TYPE " MAXIMUM - INITIAL > 100"
GO TO 900
800 CONTINUE
C ---

```

APPENDIX A  
(CONTINUED)

```

      TYPE
      DO 20 I=1,NPC
      WRITE (10,5) I
5     FORMAT (1X,"WHAT IS THE CODE FOR THE POSITIVE CONTROL #:",1X,I1,1X
1     "? (AN/6)")
      READ (11,4) (PCOM(K,I) K=1,6)
20    CONTINUE
      TYPE
      TYPE
C ---
C --- TEST COMPOUND(S)
600   CONTINUE
      ACCEPT " HOW MANY TEST COMPOUNDS ? (MAX.=5) ",NTC
      IF (NTC.GT.5) GO TO 600
      DO 30 I=1,NTC
      TYPE
      WRITE (10,6) I
6     FORMAT(1X,"WHAT IS THE CODE FOR TEST COMPOUND #:",1X,I1,1X,
1     "? (AN/6)")
      READ (11,4) (TCOM(K,I) K=1,6)

C --- GET THE REPETITION NUMBER
123   CONTINUE
      ACCEPT " WHAT IS THE REPETITION NUMBER FOR THIS COMPOUNT ? (I2) ",
1     NREP(1)
      IF (NREP(1).GT.99) GO TO 123

C ---
30    CONTINUE
      TYPE
      TYPE
C ---
C *** INPUT BROOD DATES
      TYPE
      TYPE
      TYPE "          INPUT BROOD DATES"
      DO 40 B =1,4
      WRITE (10,7) B
7     FORMAT (1X,"DATE FOR BROOD:",1X,I1,1X,"? (XXMONXX)")
      READ (11,8) (IDATE (K,B) K=1,7)
8     FORMAT (7A1)
40    CONTINUE

C ---
C *** MESSAGE
      TYPE
      TYPE "          WAIT FOR OUTPUT . . . . . "
C ---

```

APPENDIX A  
(CONTINUED)

```

C *** MAKE LABELS FOR EACH BROOD
C --- OPEN FILE 'DLABELS'
      DELETE "DLABELS"
      OPEN 1, "DLABELS", ATT = "SOP"
      DELETE "BCARDS"
      OPEN 3, "BCARDS", ATT = "SOP"

C ---
C *** OUTPUT RUN NUMBER
      WRITE (10,9) NR
      9   FORMAT (/ ,10X,"RUN:",1X,I4,/)
C ---
C *** LABELS MADE BY BROOD
      DO 50 B = 1,4
      WRITE (10,10) B
      10  FORMAT (/ ,10X,"BROOD:",1X,I1)
C ---
C --- NEGATIVE CONTROL LABELS
      K = "C"
      DO 60 I=1,NRC
      CALL NBLABELS (GSN,NR,K,I,1A,UR,CORR(1,I),B,LDATE(1,B))
      60  CONTINUE
C ---
C --- POSITIVE CONTROL LABELS
      K = "P"
      DO 70 I=1,NPC
      CALL NBLABELS (GSN,NR,K,I,1A,UR,PCOR(1,I),B,LDATE(1,B))
      70  CONTINUE
C ---
C --- TEST COMPOUND LABELS
      K = "T"
      DO 80 I=1,NTC
      CALL NBLABELS (GSN,NR,NREP(1),K,I,1A,UR,TCON(1,I),B,LDATE(1,B))
      80  CONTINUE
C ---
      TYPE
      50  CONTINUE
C ---
      CLOSE 1
      CLOSE 3

C *** END MESSAGE
      TYPE
      TYPE " ALL LABELS AND CARDS ARE READY IN 'DLABELS' AND 'BCARDS'"
      TYPE
C ---
      END

```

# APPENDIX B

## NBLABELS.FR

```

C --- SUBROUTINE NAME:          NBLABELS.FR
C --- WRITTEN BY:              WARNER JEDERBERG, CUTANEOUS HAZARDS
C --- SUBROUTINES:             NONE
C --- INPUT:                   GSN,Nr,K,P,L,M,LCOM,N,DATE
C ---                          GSN = GLP STUDY NUMBER
C ---                          Nr = RUN NUMBER
C ---                          K = "C", "P",
C ---                          L = INITIAL VALUE
C ---                          M = MAXIMUM VALUE
C ---                          LCOM = COMPOUND CODE
C ---                          N = BROOD #
C ---                          P = SUBGROUP
C ---                          IDATE = PROPER DATE PER BROOD
C --- FORMAT FOR INPUT:       SEE SLRLABELS.FR

C --- PURPOSE:                 GENERATES LABELS AND CARDS FOR THE NEGATIVE AND
C --- POSITIVE CONTROLS FOR THE SLRL DROSOPHILA ASSAY. SUPPORTS THE
C --- MAIN PROGRAM SLRLABELS.FR ALONG WITH NTBLABELS.FR.
C ---
SUBROUTINE NBLABELS (GSN,Nr,K,P,L,M,LCOM,N,DATE)
INTEGER I,P
DIMENSION IDATE(7), LCOM(6)
DO 10 I=L,M
WRITE (1,3) GSN,Nr,K,P,I,N,DATE,LCOM

3  FORMAT (1X,"GLP STUDY NO.",1X,I5,/,
11X,"RUN:",1X,I4,/,1X,A1,I1,"-",14,2X,"BR:",
21X,I1,2X,7A1,/,
21X,"COMPOUND CODE:",2X,6A1,/,1X,"NOTES:",/)
C --- PREPARE CORRESPONDING CARD
WRITE (3,7) GSN,LCOM,K,P,I,N,DATE,Nr

7  FORMAT (1X,"GLP STUDY NO.",1X,I5,/,
116X,"COMPOUND CODE:",2X,6A1,/,
22X,A1,I1,"-",14,2X,"BR:",1X,I1,2X,7A1,11X,"RUN:",1X,I4,/,
32X,"-----",/,
42X,"F2 CROSS          MEDIUM BATCH #:",/,
52X,"DATE:          INITIALS",/,
62X,"FAILURES      LETHALS      NONLETHALS",/,
72X,"-----",/,
82X,"F3 CROSS          MEDIUM BATCH #:",/,
92X,"DATE:          INITIALS",/,
12X,"FAILURES      LETHALS      NONLETHALS",/,
22X,"-----",/

```

APPENDIX B  
(CONTINUED)

```
52X,"NOTED:",//)
10 CONTINUE
C --- STORE LAST SEQUENCE # IN "LDLABEL"
    LL = 1-1
    OPEN 2, "LDLABEL", ATT = "SO"
    WRITE (2,5) NR,LL
5   FORMAT (2(1X,I4))
    CLOSE 2
C --- USER NOTE
    WRITE (10,6) K,P,LCOM
6   FORMAT (/ ,1X,"LABELS AND CARDS FOR",1X,A1,I1,"-",
16A1,1X,"ARE IN 'DLABELS' AND 'BCARDS'")
    RETURN
END
```

# APPENDIX C

## NTBLABELS.FR

```

C --- SUBROUTINE NAME:      NTELABELS.FR
C --- WRITTEN BY:          WARREN JEDERBERG, CUTANEOUS HAZARDS
C --- SUBROUTINES:         NONE
C --- INPUT:               GSN,Nr,NREP,K,P,L,X,LCOM,N,DATE
C ---                     GSN = GLP STUDY NUMBER
C ---                     Nr = RUN NUMBER
C ---                     NREP = REPETITION NUMBER

C ---                     K= "T"
C ---                     L = INITIAL VALUE
C ---                     M = MAXIMUM VALUE
C ---                     LCOM = COMPOUND CODE
C ---                     N = CROSS #
C ---                     P = SUBGROUP
C ---                     DATE = PROPER DATE PER CROSS
C --- FORMAT FOR INPUT:    SEE SRLABELS.FR

C --- PURPOSE:             SUBROUTINE TO GENERATE LABELS AND CARDS FOR THE
C --- TEST COMPOUNDS USED IN THE SRL DROSOPHILA ASSAY. SUPPORTS THE
C --- PROGRAM SRLABELS.FR ALONG WITH NPLABELS.FR.
C ---
SUBROUTINE NTELABELS (GSN,Nr,NREP,K,P,L,X,LCOM,N,DATE)
INTEGER I,P
DIMENSION DATE(7), LCOM(6)
DO 10 I=L,M
WRITE (1,3) GSN,Nr,NREP,K,P,I,N,DATE,LCOM

3   FORMAT (1X,"GLP STUDY NO.",1X,15,/,
11X,"RUN:",1X,14,5X,"REP. #:",1X,12,/,1X,A1,11,"-",14,2X,"FA:",
11X,11,2X,7A1,/,
21X,"COMPOUND CODE:",2X,6A1,/,1X,"NOTES:",/)
C --- PREPARE CORRESPONDING CARD
WRITE (3,7) GSN,NREP,LCOM,K,P,I,N,DATE,Nr

7   FORMAT (1X,"GLP STUDY NO.",1X,15,4X,"REP. #:",1X,12,/,
118X,"COMPOUND CODE:",2X,6A1,/,
22X,A1,11,"-",14,2X,"Ba:",1X,11,2X,7A1,11X,"RUN:",1X,14,/,
52X,"-----",/,
42X,"F2 CROSS          MEDIUM BATCH #:",/,
52X,"DATE:          INITIALS",/,
62X,"FAILURES      LETHALS      NON-LETHALS",/,
72X,"-----",/,
82X,"F3 CROSS          MEDIUM BATCH #:",/

```

APPENDIX C  
(CONTINUED)

```

9XX,"DATE: _____ INITIALS _____",//,
1XX,"FATALITIES _____ LETHALS _____ NONLETHALS _____",/,
2XX,"-----",/,
3XX,"NOTES:",//)
10 CONTINUE
C --- STORE LAST SEQUENCE # IN "LDLABEL"
  LL = 1-1
  OPEN 2, "LDLABEL", ATT = "SO"
  WRITE (2,5) NR,LL
  FORMAT (2(1X,I4))
  CLOSE 2
C --- USER NOTE
  WRITE (10,6) K,P,LCOM
  6 FORMAT (/,1X,"LABELS AND CARDS FOR",1X,A1,I1,"-",
  10,1X,"ARE IN 'DLABELS' AND 'BCARDS'")
  RETURN
END

```

END

DATE  
FILMED

6-82

DTIC